

# Cell Scrapers



"Scraper" Blade in use



"Lifter" Blade in use

- Cell Scrapers are designed to gently harvest cells from tissue culture flasks, dishes or bottles
  - Soft, pliable blade won't damage cells or scratch surfaces
- 1 "Scraper" Blade is positioned parallel to the handle for windshield wiper like scraping (PN 229310 & 229315)
  - 2 "Lifter" Blade is positioned perpendicular to the handle for use with a raking motion (PN 229311 & 229316)
  - 3 Pivoting Blade Cell Scraper features free rotating blade that adjusts to the desired position, providing access to all corners of the flask or tissue culture vessel (PN 229312)

## CELL SCRAPERS

S, R/D F, NP, USP, ABS/TPE

Part No.	Product Description	Handle Length	Blade Length	Blade Orientation	Pack Description	Qty/Case
229310*	Cell Scraper	25cm	20mm	Scraper (parallel)	Individually Wrapped	100
229311	Cell Scraper	25cm	20mm	Lifter (perpendicular)	Individually Wrapped	100
229312	Cell Scraper	30cm	20mm	Pivoting	Individually Wrapped	100
229315	Cell Scraper	39cm	30mm	Scraper (parallel)	Individually Wrapped	100
229316	Cell Scraper	39cm	30mm	Lifter (perpendicular)	Individually Wrapped	100

\*Most popular option

877.345.4620 • [www.celltreat.com](http://www.celltreat.com) 27